

**JAVALYTICS PRODUCTS
COMPARISON**

	JAV-RDA-U <i>"Up Looker"</i>	JAV-RDA-D <i>"Down Looker"</i>	JAV-RDA-DN <i>"Down Looker"</i>
US List Price	\$7,600 USD	\$10,400 USD	\$9,401 USD
Network (JAV-NET)	*Optional	Included	*Optional
Glassware	Must maintain petri dishes: keep clean and unbroken. Polystyrene petri dishes are available, but not green. Some locations can not use glass.	Not utilized. The samples are prepared in a solid plastic tray.	Not utilized. The samples are prepared in a solid plastic tray.
Ground Coffee	Samples self-level on bottom of petri dish. Careful attention for ground coffee, as settling of fines will move toward lighter reading. Reproducibility enhanced with repeatable grinding and settling technique.	Good technique important for reproducibility. Sample screeding technique similar to smoothing concrete. Again, grind size control is important. Optional sample loader/leveler available.	Good technique important for reproducibility. Sample screeding technique similar to smoothing concrete. Again, grind size control is important. Optional sample loader/leveler available.
Whole Beans	Routine, as beans settle flat to dish. May want to average more than one cup if obvious color or chaff variations are visible.	A bit more difficult to reproducibly get "level". May want to average more than one cup if obvious color or chaff variations are visible	A bit more difficult to reproducibly get "level". May want to average more than one cup if obvious color or chaff variations are visible.

JAVALYTICS instruments are very stable and offer great reproducibility. Variation is primarily due to sample preparation.

** The JAV-NET option is \$999.*

